

## National Guidance for Oesophageal Doppler Monitoring (ODM)

### Background

The oesophageal Doppler monitor (ODM) was first released onto the market place in 1992 with the first Randomised Controlled Trial (RCT) being published in 1995.

Since its release the ODM has been used on over 1,000,000 patients in over 40 countries around the world.

ODM-guided fluid management has repeatedly shown reductions in post-op complications and length of hospital stay across 24 published RCT's in multiple different types of surgery.

These trials have been subject to a number of meta-analyses and systematic reviews:

- A health technology assessment (HTA) published by ECRI for the Agency for Healthcare Research & Quality (AHRQ);
- NICE recommendation (MTG3) for over 800,000 medium and high-risk surgery patients;
- Included within NICE pathway for 'Major haemorrhaging in hospital';
- Multiple meta-analyses re ODM in abdominal surgery published;
- Multiple meta-analysis re ODM across all surgical applications published;
- Multiple meta-analysis re ODM and GDTF in all hospital settings;
- A report by the Centers for Medicare and Medicaid Services based upon the ECRI HTA;
- A report by the UK Centre for Evidence-based Purchasing (CEP);
- An HTA by the Aberdeen University Health Economic Research Unit commissioned by the NHS National Institute for Healthcare Research (NIHR);
- Assessment of ODM evidence as Level 1A for abdominal surgery and Level 1B for orthopaedic surgery for UK consensus fluid guidelines;
- 24 Randomised Clinical Trials across different surgical applications (colorectal, urology, gynecology, orthopaedic, hepatic) in adult and paediatric patients
- Over 1,000 peer reviewed publications of utilising the ODM in the operating room, intensive care and emergency room.

All the above have concluded positively on the benefits of ODM. The AHRQ/ECRI report concluded that, 'for a large group of surgical patients, no further clinical trials on ODM were necessary'.

### Main Recommendations

#### UK

National Institute For Health and Care Excellence (NICE)

<https://www.nice.org.uk/guidance/MTG3/chapter/1-recommendations>

The case for adopting the CardioQ-ODM in the NHS, when used as described in 1.2, is supported by the evidence. There is a reduction in post-operative complications, use of central venous catheters and in-hospital stay (with no increase in the rate of re-admission or repeat surgery) compared with conventional clinical assessment with or without invasive cardiovascular monitoring. The cost saving per patient, when the CardioQ-ODM is used instead of a central venous

catheter in the peri-operative period, is about £1100 based on a 7.5-day hospital stay.

1.2 The CardioQ-ODM should be considered for use in patients undergoing major or high-risk surgery or other surgical patients in whom a clinician would consider using invasive cardiovascular monitoring.

### **NICE Pathway Major haemorrhaging in hospital'**

<https://pathways.nice.org.uk/pathways/trauma>

### **USA**

#### **Medicare National Coverage Decision (NCD)**

<https://www.cms.gov/outreach-and-education/medicare-learning-network-mln/mlnmattersarticles/downloads/MM8330.pdf>

The Centres for Medicare and Medicaid Services (CMS) provides coverage for use of Esophageal Doppler Monitor (EDM). CMS determined that EDM is “reasonable and necessary” for all Medicare beneficiaries who meet the criteria for the following patient populations:

- Monitoring of cardiac output for ventilated patients in the ICU; and
- Operative patients with a need for intra-operative fluid optimization

The 2015 national average Medicare payment level for G9157 is \$97.97 (based on 2.74 RVUs). Professional services associated with EDM are separately payable under the Medicare Physician Fee Schedule (MPFS). Medicare adjusts professional payment based on regional differences in practice costs, thus payment may vary based on physician location. Other payers, including private payers, may base their payments on Medicare payment rates.

#### **AHRQ Technology Assessment**

[https://www.ecri.org/Resources/EPC\\_Sample\\_Reports/Esophageal Doppler Ultrasound Based Cardiac Output Monitoring.pdf](https://www.ecri.org/Resources/EPC_Sample_Reports/Esophageal_Doppler_Ultrasound_Based_Cardiac_Output_Monitoring.pdf)

Esophageal Doppler Ultrasound-Based Cardiac Output Monitoring for Real-Time Therapeutic Management of Hospitalized Patients

### **Spain**

**Spanish Ministry of Health, Social Services and Equality**

[http://portal.guiasalud.es/contenidos/iframes/documentos/opbe/2015-07/ViaClinica-RICA\\_English.pdf](http://portal.guiasalud.es/contenidos/iframes/documentos/opbe/2015-07/ViaClinica-RICA_English.pdf) (English Version)

National guidance for Spain to conduct all surgery within enhanced recovery protocols. Oesophageal Doppler is the recommended technology to use within the protocol.

#### Page 51 – Summary and Recommendations

60. The use of SV (Stroke Volume) or of SVV (Stroke Volume Variation) by monitoring is recommended to guide the intraoperative administration of fluids.  
Strong recommendation +. High Level of Evidence.
62. Restrictive continuous fluid perfusion must be maintained in order to avoid fluid overload.  
Strong recommendation +. High Level of Evidence.
63. Intraoperative hypotension must be treated with vasopressors.  
Strong recommendation +. High Level of Evidence
64. A mean blood pressure range of 70 mmHg must be established  
Strong recommendation +. Moderate Level of Evidence
65. A CI of > 2.5 l/min/m<sup>2</sup> must be maintained, using inotropes in cases where there is no response to volume.  
Strong recommendation +. Moderate Level of Evidence.
66. **Monitoring with oesophageal Doppler is preferred**, or else methods based on validated pulse contour analyses.  
Strong recommendation +. Moderate Level of Evidence.

#### France

##### Société Française d'Anesthésie et de Réanimation ('SFAR')

<http://sfar.org/strategie-du-remplissage-vasculaire-perioperatoire-2/>

##### (National guidelines for peri-operative optimisation during surgery)

The guidelines are divided into 15 recommendations regarding fluid management, with 3 directly relevant to the use of the ODM

Each recommendation made within the document is graded in a binary system:

- High: You must do or not do (**GRADE1+** or **GRADE1-**)
- Low: It is possible to do or not do (**GRADE2+** or **GRADE2-**)

All 3 recommendations relating directly to ODM fluid management received a GRADE1+ rating.

This means that they must all be carried out in the case of high-risk patients, as the evidence level is high. It has been assessed by the SFAR that future evidence is unlikely to change any conclusions drawn from current evidence.

The three relevant recommendations provided are as follows:

**Q1:** Filling guided by the measurement of the systolic ejection volume of (SV) reduces post-operative morbidity and duration of stay.

**A1:** Fluid titration for high-risk surgical patients should be guided by stroke volume to reduce post-operative morbidity, get earlier return to oral feeding and reduce length of stay. **GRADE1+.**

**Q2:** Should we interrupt the filling in the absence of the SV increase?

**A2:** It is recommended to discontinue filling if SV does not increase. **GRADE1+**

**Q3:** Should we regularly reassess the SV?

**A3:** Reassess SV regularly and the SV response to a fluid challenge, especially during periods of haemodynamic instability. **GRADE1+.**

## Poland

Polskie wytyczne okołoperacyjnego leczenia płynami

<https://ojs.kardiologiapolska.pl/kp/article/download/KP.2014.0193/7958>.

Page 8 - National recommendations for the use of SV Optimisation using the ODM

## Enhanced Recovery After Surgery (ERAS)

Development and Feasibility Study of an Algorithm for Intraoperative Goal Directed Haemodynamic Management in Noncardiac Surgery

A Feldheiser, P Conroy, T Bonomo, B Cox, T Ruiz Garces and C Spies on Behalf of the Anaesthesia Working Group of the Enhanced Recovery After Surgery (ERAS®) Society

'Preload optimization guided by stroke volume, as measured by oesophageal Doppler ultrasonography, was almost always associated with a significant reduction in length of hospital stay and postoperative morbidity (Table 2). In contrast to these studies, investigations based on the measurement of stroke volume by the Vigileo™ monitor/FloTrac™ sensor (Edwards Lifesciences, Irvine, CA, USA) or by the lithium indicator dilution cardiac output (LiDCO™; Lidco, Cambridge, UK) system did not uniformly improve outcome (Table 2).

'There is no evidence that the use of a pulse contour method to guide intraoperative fluid therapy by measuring stroke volume is equivalent to guiding fluid administration by oesophageal Doppler ultrasonography'.

ERAS Recommendation

'Oesophageal Doppler ultrasonography was chosen as the preferred method of monitoring intraoperative stroke volume due to the broader evidence base in this context.'

## Enhanced Recovery Partnership Guidance UK

<http://perioperativemedicinejournal.biomedcentral.com/articles/10.1186/2047-0525-1-2>

'The Enhanced Recovery Partnership fully supports the use of intra-operative fluid management technologies to deliver individualised goal directed fluid therapy. This is recommended in the 2012-13 NHS Operating Framework, in the Innovation, Health and Wealth Review and in NICE Guideline MTG3.'

## Oesophageal Doppler Monitoring (ODM)

### Published Randomised Controlled Trials (RCT's):-

1. Mythen MG, Webb AR. Perioperative plasma volume expansion reduces the incidence of gut mucosal hypoperfusion during cardiac surgery. *Arch Surg* 1995; 130: 423-9
2. Sinclair S, James S, Singer M. Intraoperative intravascular volume optimisation and length of hospital stay after repair of proximal femoral fracture: randomised controlled trial. *Br Med J* 1997; 315: 909-12
3. Conway DH, Mayall R, Abdul-Latif MS, Gilligan S, Tackaberry C. Randomised controlled trial investigating the influence of intravenous fluid titration using oesophageal Doppler monitoring during bowel surgery. *Anaesth* 2002; 57: 845-9
4. Gan TJ, Soppitt A, Maroof M, El-Moalem H, Robertson KM, Moretti E, Dwane P, Glass PSA. Goal-directed intraoperative fluid administration reduces length of hospital stay after major surgery. *Anesthesiol* 2002; 97: 820-6
5. Venn R, Steele A, Richardson P, Poloniecki J, Grounds M, Newman P. Randomized controlled trial to investigate influence off the fluid challenge on duration of hospital stay and perioperative morbidity in patients with hip fractures. *Br J Anaesth* 2002; 88(1): 65-71
6. McKendry, M., McGloin, H., Saberi, D., Caudwell, L., Brady, A.R., and Singer, M. *Randomised controlled trial assessing the impact of a nurse delivered, flow monitored protocol for optimisation of circulatory status after cardiac surgery. BMJ, 2004. 329(7460): p. 258.*
7. Wakeling HG, McFall MR, Jenkins CS, Woods WGA, Miles WFA, Barclay GR, Fleming SC. Intraoperative oesophageal Doppler guided fluid management shortens postoperative hospital stay after major bowel surgery. *Br J Anaesth* 2005; 95(5): 634-642
8. Noblett SE, Snowden CP, Shenton BK, Horgan AF. Randomized clinical trial assessing the effect of Doppler-optimized fluid management on outcome after elective colorectal resection. *Br J Surg* 2006; 93: 1069-76
9. Chytra, I., Pradl, R., Bosman, R., Pelnar, P., Kasal, E., Zidkova, A., *Esophageal Doppler-guided fluid management decreases blood lactate levels in multiple-trauma patients: a randomized controlled trial. Crit Care, 2007. 11(1): p. R24*

10. Senagore AJ, Emery T, Luchtefeld M, Kim D, Dujovny N, Hoedema R. Fluid management for laparoscopic colectomy: a prospective, randomized assessment of goal-directed administration of balanced salt solution or hetastarch coupled with an enhanced recovery program. *Dis Colon Rectum* 2009; 52(12): 1935-40
11. Challand C, Struthers R, Sneyd JR, Erasmus PD, Mellor N, Hosie B, Minto G. Randomized controlled trial of intraoperative goal-directed fluid therapy in aerobically fit and unfit patients having major colorectal surgery. *Br J Anaesth* 2012; 108(1): 53-62
12. Pillai P, McEleavy I, Gaughan M, Snowden C, Nesbitt I, Durkan G, Johnson M, Cosgrove J, Thorpe. A double-blind randomized controlled clinical trial to assess the effect of Doppler optimized intraoperative fluid Management on outcome following radical cystectomy. *J Urol* 2012; 186(6): 2201-6
13. Brandstrup B, Svendsen PE, Rasumssen M, Belhage B, Rodt SÅ, Hansen B, Møller DR, Lundbech LB, Andersen N, Berg V, Thomassen N, Andersen ST, Simonsen L. Which goal for fluid therapy during colorectal surgery is followed by the best outcome: near maximal stroke volume or zero fluid balance? *Br J Anaesth* 2012; 109(2): 191-9  
12
14. Srinivasa S, Taylor MGH, Singh PP, Yu T-C, Soop M, Hill AG. Randomized clinical trial of goal-directed fluid therapy within an enhanced recovery protocol for elective colectomy. *Br J Surg* 2013; 100(1): 66-74
15. Zakhaleva J, Tam J, Denoya PI, Bishawi M, Bergamaschi R. The impact of intravenous fluid administration on complication rates in bowel surgery within an enhanced recovery protocol: a randomized controlled trial. *Colorectal Dis* 2013; 15: 892-9
16. McKenny M, Conroy P, Wong A, Farren M, Gleeson N, Walsh C, O'Malley C, Dowd N. A randomised prospective trial of intra-operative oesophageal Doppler-guided fluid administration in major gynaecological surgery. *Anaesthesia* 2013
17. El Sharkawy OA, Refaat EK, Ibraheem AEM, Mahdy WR, Fayed NA, Mourad WS, Abd Elhafez HS, Yassen KA. Transoesophageal Doppler compared to central venous pressure for perioperative hemodynamic monitoring and fluid guidance in liver resection. *Saudi J Anaesth* 2013; 7(4): 378-86
18. Phan TD, D'Souza B, Rattray MJ, Johnston MJ, Cowie BS. A randomised controlled trial of fluid restriction compared to oesophageal Doppler guided goal directed therapy in elective major colorectal surgery within an Enhanced Recovery After Surgery program. *Anaesth Intensive Care* 2014 Nov; 42 (6);752-60
19. Eman Sayed Imbrahim, Taha Aid Yassein, Wesam Saber Morad. The beneficial values of transoesophageal Doppler in intraoperative fluid guidance versus standard clinical monitoring parameters in infants undergoing Kasai operation. *Middle East Journal Anaesthesiology* 2015 Jun; 23(2); 205-11
20. Picard J, Bedague D, Bouzat P, Ollinet C, Albaladejo P, Bosson JL, Payen JF. Oesophageal Doppler to optimize intraoperative haemodynamics during prone position. A randomized controlled trial. *Anaesth Crit Care Pain Med*. 2016 Aug;35(4):255-60. doi: 10.1016/j.accpm.2015.12.011.

21. Kaufmann KB, Steil L, Ulbrich F, Kaifi JT, Hauschke D, Loop T, Goebel U. Oesophageal Doppler guided goal-directed haemodynamic therapy in thoracic surgery - a single centre randomized parallel-arm trial. *Br J Anaesth.* 2017 Jun 1;118(6):852-861. doi: 10.1093/bja/aew447.
22. Feldheiser A, Pavlova V, Weinmann K, Hunisicker O, Stockmann M, Koch M, Giebels A, Wenecke KD, Spies CD. Haemodynamic Optimization by Oesophageal Doppler and Pulse Power Wave Analysis in Liver Surgery: A Randomised Controlled Trial. *PLoS One.* 2015 Jul 17;10(7):e0132715 doi: 10.1371/journal.pone.0132715
23. Calvo-Vecino JM, Ripollés-Melchor J, Mythen MG, Casans-Francés R, Balik A, Artacho JP, Martínez-Hurtado E, Serrano Romero A, Fernández Pérez C, Asuero de Lis S; FEDORA Trial Investigators Group. Effect of goal-directed haemodynamic therapy on postoperative complications in low-moderate risk surgical patients: a multicentre randomised controlled trial (FEDORA trial). *Br J Anaesth.* 2018 Apr;120(4):734-744.
24. Szturz P, Folwarczny P, Kula R, Neiser J, Ševčík P, Benes  
Multi-parametric functional hemodynamic optimization improves postsurgical outcome after intermediate risk open gastrointestinal surgery, a randomized controlled trial. *Minerva Anesthesiol.* 2018 May 11. doi: 10.23736/S0375-9393.18.12467-9
25. Muhlbacher J, Luf F, Zotti O, Herkner H, Fleischmann E, Kabon B  
Effect of Intraoperative Goal-Directed Fluid Management on Tissue Oxygen Tension in Obese Patients: a Randomized Controlled Trial. *Obesity Surgery* 2020 November 27 1129–1138(2021)
26. Respiratory and Hemodynamic Effects of Prophylactic Alveolar Recruitment During Liver Transplant: A Randomized Control Trial. *Experimental and Clinical Transplantation* May 2021

Last Updated 1<sup>st</sup> June 2021