(UK) - Instructions For Use (IFU): Kinder Oesophageal Doppler Probe

1. Intended Use:

The Kinder Oesophageal Doppler Probe is used with the CardioQP™/CardioQ ODM™/CardioQP ODM™ Monitoring System.

Clinical benefit- Proven effective for haemodynamic optimisation to reduce complications and shorten length of hospital stay

Indications- The Deltex Medical CardioQ-ODM+ is indicated for use as a fluid management and cardiac output monitoring system

The CardioQP/CardioQ ODM/ CardioQP ODM contains a Nomogram which can be used to calculate the cardiac output in paediatric patients.

The paediatric nomogram is valid for patients -

Age 0 to 15 yrs.

Weight 3 kg to 60 kg (33 lb to 132 lb).

Height 50 cm to 170 cm (20" to 67").

The Kinder Oesophageal Doppler Probe is for oral use in patients under general anaesthesia.

Before use read these instructions in conjunction with the Operating Handbook.

2. Warnings - Do Not Use in the Following Situations - Patient Harm May Result:

- Do not use after the expiry date printed on the shelf box and unit pack labels.
- Do not use undue force upon insertion. Remove if difficulty arises during insertion and seek advice.
- Use only by trained paediatricians or under such supervision.
- Not recommended for use in patients under going intra-aortic balloon pumping.
- The probe may be unreliable during periods of aortic cross clamping.
- Consider coagulation status of patient prior to insertion of probe.
- Down's syndrome children may have intrinsic narrowing of their hypopharyngeal structures preventing insertion of the probe.
- It is not recommended that the KDPn probes are left in place for extended periods of time.
- Any serious incident that occurs while using this product shall be reported to Deltex Medical Ltd and your member state competent authority.
- MR unsafe, do not use this equipment in the MRI scan room

Single patient use only. Risk of cross contamination if probe is used on more than one patient

3. Cautions:

Do not use if the pouch has been unintentionally opened, or if the probe or pouch show signs of damage. The probe shall be disposed of according to local procedures

Exercise caution when handling. The probe contains an internal spring that will cause the probe to uncoil when released from the Pouch.

Open the pouch at the chevron seal end.

Do not use scissors or sharp tools to open the pouch.

Consider insertion under direct vision.

Clinicians are advised to take into consideration the potential of increased risks of pressure on the oropharynx when used in combination with other equipment.

Clinicians should consider the removal of the probe between periods of measurement.

Should the probe be removed and stored for same patient reuse, refer to local guidelines for cleaning and storage.

4. Contraindications:

Do not use on patients weighing less than 3 kg.

Do not insert nasally.

Do not use with carcinoma of the pharynx; larynx; oesophagus.

Do not use with aneurysms of the thoracic aorta.

Do not use with proximal coarctation of the aorta.

Do not use with tissue necrosis of the oesophagus.

Do not use in close proximity to laser surgery

Not for use in patients with pharyngo-oesophago-gastric pathology and/or severe bleeding diatheses

5. General Information:

The Doppler Probes are supplied STERILE.

The Probe can only be used with a Deltex Medical CardioQP/CardioQ ODM/ CardioQP ODM Monitoring System.

Probe output specifications is documented within the ODM Monitor handbook

6. Removal of the Probe from the Unit Pack:

Open the pouch at the chevron seal end and, holding the Probe and, withdraw it from the pouch.

Allow the Probe to uncoil, keeping it free from contact with any object that may compromise its state of cleanliness.

7. Insertion of the Probe into the Patient:

Insertion should be carried out by an appropriately trained and qualified paediatrician or by a person under such supervision.

The Probe must be inserted orally. Note the Warnings above.

The Probe connector should be plugged into the Patient Interface Cable connector and the correct nomogram details entered.

Apply a water-based lubricant to the tip of the Probe. This aids insertion and signal acquisition.

Insert the Probe with the bevelled edge of the tip facing the back of the throat so that when in the oesophagus, the bevelled edge faces backwards and downwards towards the descending aorta.

Excessive force should not be used when inserting the Probe. Remove if undue difficulty is experienced.

Warning - Risk of patient harm if forced.

The numbered depth markers on the Probe shaft are for use as a guide for oral placement. Signal acquisition and depth of insertion will depend upon patient height. See Table 1 as a guide to expected insertion depths; use both direct visual and audible focusing to achieve the optimum ultrasound signal.

Refer to the Operating Handbook for positioning the probe and operating the Monitoring System.

The total probe usage time remaining is displayed when the probe is connected.

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Proizvodi Deltex Medical Ltd, Terminus Road, Chichester, PO19 8TX, UK
Consultancy Services Ltd.
medReg GmbH